# The Reduction in Water Consumption of Sterilizer Equipment Resulting From the Installation of Water-Mizer<sup>™</sup> Systems

A Technical Evaluation Prepared by:

TDK Consulting Services 12921 NW 78<sup>th</sup> Terrace Parkville, MO 64152





I certify that the information herein has been prepared under my direct supervision and the contents are true and correct.

## 1 Executive Summary

TDK Consulting Services was contracted by Continental Equipment Company Inc. to evaluate the reduction in water consumption of sterilizer equipment resulting from the installation of Water-Mizer<sup>TM</sup> systems. The Water-Mizer is a tempering device that mixes cold water with hot water discharged from sterilization equipment to reduce the discharged water temperature. There is an additional system available that can be added to the original Water-Mizer that captures and re-uses water that is necessary to produce the vacuum for this equipment. Both of these systems were evaluated.

The testing consisted of installing Water-Mizer equipment on two of three sterilizers of the same make and model, leaving the third sterilizer as a control, and observing the operation of all three units during comparable operation. The testing configurations are shown in Table 1-1. Four Sterilization Cycles were observed for each unit during a one-day period. Water consumption and other observations were recorded.

| Installations |   |  |  |
|---------------|---|--|--|
| Unit 1        | Water-Mizer with Recirculation System Installed |  |  |
| Unit 2        | Water-Mizer Installed                           |  |  |
| Unit 3        | No Water-Mizer Equipment Installed              |  |  |

| Table 1-1 resund Configurations | Table 1-1 | Testing | Configurat | ions |
|---------------------------------|-----------|---------|------------|------|
|---------------------------------|-----------|---------|------------|------|

The testing was successfully completed with no operational problems. Results are summarized in Table 1-2. Calculations with report data are contained in Appendix G and are referenced in report text with the use of footnotes.

| Table 1-2 Test Results, gallons/cycle | Э |
|---------------------------------------|---|
|---------------------------------------|---|

|         |        | , 0    |        |
|---------|--------|--------|--------|
|         | Unit 1 | Unit 2 | Unit 3 |
| Test 1  | 62     | 253    | 356    |
| Test 2  | 74     | 281    | 364    |
| Test 3  | 86     | 287    | 372    |
| Test 4  | 80     | 289    | 418    |
| Average | 76     | 278    | 378    |



As shown in Table 1-2, both units equipped with Water-Mizer equipment used less water than the sterilizer with no Water-Mizer equipment installed.

Conclusions from this testing are as follows:

- The installation of the Water-Mizer including the recirculation system can reduce Sterilizer Cycle water consumption between 73<sup>1</sup> and 80<sup>2</sup> percent. Comparable savings are anticipated for other cycles such as the Dart Test and Leak Test.
- The installation of the Water-Mizer including the recirculation system can reduce total sterilizer water consumption by more than 80 percent. Although not directly tested in this evaluation, these reductions could exceed 90 percent<sup>3</sup>.
- The addition of the recirculation system to the Water-Mizer can reduce Sterilizer Cycle water consumption by 73<sup>1</sup> percent. Comparable savings are anticipated for other cycles such as the Dart Test and Leak Test.
- The installation of the Water-Mizer can drastically reduce standby losses.
- Depending upon the frequency of use, the manual valve settings, and the piping configuration of the condensate drains, sterilizer standby water losses will likely exceed the amount of water consumed during operational cycles (over prolonged periods).
- The cost of operating the electric pump associated with the Water-Mizer recirculation system is negligible when compared to the cost of water saved with its use.
- Considering water and sewer costs to be 5.0 \$/kgal, electricity costs to be 6.0 ¢/kwh, and the total sterilizer water consumption savings to be 90 percent, the installation of the Water-Mizer equipment on Unit 3 at this location would save the owner 9,000 \$/yr in utility costs<sup>4</sup>. This savings represents a typical application of a Water-Mizer installation.



## 2 Testing Parameters

The testing was performed at a major metropolitan hospital located in Denver, Colorado on April 5, 2005. The hospital facility staff was aware of and approved the testing exercise. On March 8, 2005, the Water-Mizer was installed on one of the hospital's sterilizers. On April 4, 2005, the Water-Mizer (including the recirculation system) was installed and commissioned on another of the hospital's sterilizers. A third sterilizer with no Water-Mizer equipment was also tested as a control device. The three sterilizers are located adjacent to one another and are of same make and model (AMSCO 3043 with Stage III Controller). See Appendix F for photographs of the units. Unit 1 was re-manufactured in 1995 and Units 2 and 3 were remanufactured in 1996. The three sterilizers involved in the testing are used on a regular basis to prepare surgical instruments and their operating schedules are considered representative for installations of this type.

This testing was performed to evaluate water consumptions only and did not evaluate the functional performance or design integrity of the Water-Mizer systems. TDK Consultants are not contracted to consult on design issues of this product.

### **Sterilizer Cycles and Water Consumptions**

The three sterilizers are operated with three different cycles; a Dart Test, a Leak Test, and the Sterilization Cycle. The Dart Test is used to evaluate the sterilizer's ability to operate under a vacuum. The Leak Test is used to evaluate the integrity of the chamber seal. The Sterilizer Cycle is the functional cycle used to prepare surgical instruments. Printouts of operational logs for each of these three cycles are included in Appendix A. The programmed sequence of these three cycles is identical for each of the three machines. The sequences are not altered, except in rare occasions when the Sterilization Cycle (also referred to as Cycle 1) may be changed to accommodate special conditions. According to those familiar with the use of these machines, all three sterilizers run 6 to 8 Sterilization Cycles per day, Monday through Friday and one Dart Test and one Leak Test per night, seven nights per week. On the weekends, Unit 3 is generally not used for sterilization, and Units 1 and 2 run 4 to 5 Sterilzation Cycles each day. In preparation for this testing, the cycle counts and water consumption for each of the units was recorded starting March 10, 2005; these data are included in Appendix B. Prior to the installation of the Water-Mizer equipment, the three units consumed approximately the same quantity of water for each of the three cycles, as summarized in Table 2-1. Water consumption per cycle will vary slightly by unit, chamber contents, steam conditions, domestic water temperature, and other factors.



|                   | Typical Water Consumption<br>Per Cycle, gallons | # of Annual Cycles (for all<br>3 Units, Typical) |
|-------------------|---|--|
| Leak Test         | 174   | 1,095  |
| Dart Test         | 113   | 1,095  |
| Sterilizer Cycle  | 372   | 6,417  |
| Total (estimated) | 2,701,389 (Per Year)                            | 8,607  |

Table 2-1 Pre-Water-Mizer Cycle Water Consumptions

The jacket trap drain of Unit 3 is piped directly to the sewer system (floor drain). As a result, a constant stream of water, referred to as tempering water herein, is introduced into this drain to ensure the sewer discharge water temperature does not exceed 140°F. During testing this tempering water flow rate was measured at 3.1 gpm. As a result, Unit 3 experiences stand-by losses as shown in Table 2-2. Units 1 and 2 do not experience standby losses with the Water-Mizer systems installed and the associated tempering water needle valves on these units were removed during Water-Mizer installation. Prior to Water-Mizer installations, Units 1 and 2 experienced standby losses of 2.7 and 2.8 gpm respectively.

|                   | Typical Water Consumption,<br>gpm | # of Standby Minutes<br>(Annual, Typical) |
|-------------------|-----------------------------------|---|
| Standby           | 3.1                               | 400,285                                   |
| Total (estimated) | 1,240,884 (Gallons Per Year)      |   |

#### Table 2-2 Stand-by Water Consumption of Unit 3

### **Unit Configurations**

#### Unit 1

Unit 1 was equipped with the Water-Mizer system including the recirculation system. The schematic of this testing configuration is shown in Appendix C. The recirculation system allows the water used to create the chamber vacuum, via the Bernoulli principle, to be captured and reused. The water will be recirculated with a pump as long as the temperature of the water to the tank remains below setpoint. For this test, the setpoint was  $85^{\circ}$ F. When the temperature of the water of the water to the tank rises above the setpoint, the tank inlet valve closes and the valve to the Water-Mizer tank opens, sending the water to the drain. The Water-Mizer introduces tempering water as necessary to ensure that the discharge water temperature does not exceed  $140^{\circ}$ F.

#### Unit 2

Unit 2 was equipped with the Water-Mizer without the recirculation system. The schematic of this testing configuration is shown in Appendix C. The Water-Mizer will measure the sterilizer



discharge water temperature and temper that discharge with domestic cold water as necessary to maintain a maximum outlet temperature of 140°F.

#### Unit 3

Unit 3 did not have Water-Mizer equipment installed. As discussed under *Water Consumptions* above, the difference in piping between Unit 3 and Units 1 and 2 is the disposal of the jacket condensate to drain. On Units 1 and 2 the jacket condensate is directed to the boiler plant condensate system. On Unit 3 the jacket condensate is sent to the floor drain with a constant stream of tempering water added to ensure the sewer system is not subjected to water temperatures above 140°F. As a result, the sterilizer is consuming (and disposing of) the tempering water flow rate regardless of whether the sterilizer is being used.

### **Testing Equipment**

The water consumption of each unit was measured with a positive displacement meter manufactured by Hayes Fluid Controls. The meters were installed in the domestic water supply line to each unit in anticipation of this testing. The serial numbers for the three meters were 31625014, 31625061, and 31625013, for Units 1,2, and 3 respectively. See Appendix F for photographs. The meters were manually read to the nearest gallon immediately before and after each testing cycle. Water consumption test data are presented in Section 3, Test Results.

The power consumption of the recirculation pump motor associated with the Water-Mizer installed on Unit 1 was manually read with a Fluke T5-600 electrical tester. The amperage was measured through one of the conductors with an OpenJaw current measurement device and the voltage was measured with probes inserted into the 120 volt receptacle. The amperage draw remained constant throughout testing: 6.9 amps when the recirculation system was delivering water to the drain and 7.2 amps when the recirculation system was delivering water to the recirculation tank. The voltage also remained constant throughout the test at 107 volts AC (single phase).

The temperature of the domestic water supplied to the units was measured at 49°F with a type K thermocouple. The pressure of the water was measured at 60 psig with a pressure gage installed downstream of the water consumption meter.

The conditions of the steam supply to the sterilizer and the discharge to drain were not measured.



### Schedule

Four Sterilization Cycles were tested on each unit. The schedule of testing is provided in Table 2-3. Dart Test and Leak Test cycles were not evaluated in this effort.

|        | Unit 1      | Unit 2      | Unit 3      |
|--------|-------------|-------------|-------------|
| Test 1 | 09:19-10:04 | 09:08-09:53 | 08:51-09:36 |
| Test 2 | 10:35-11:31 | 10:16-11:05 | 09:45-10:31 |
| Test 3 | 12:33-13:30 | 13:36-14:26 | 10:56-11:43 |
| Test 4 | 15:06-16:09 | 15:40-16:30 | 14:51-15:43 |

#### Table 2-3 Testing Schedule -- April 5, 2005

Test 1 on each unit was run with the chamber empty. Tests 2 - 4 on each unit were run with the chamber loaded by hospital staff. The operations of the sterilizers during testing are documented in Appendix D, Test Log Printouts.



## 3 <u>Test Results</u>

The testing was accomplished without incident. All three units successfully completed four Sterilization Cycles without operational difficulty. The water consumption of each test is provided in Table 3-1. Additional graphical representations of test data are presented in Appendix E.

|         | Unit 1 | Unit 2 | Unit 3 |
|---------|--------|--------|--------|
| Test 1  | 62     | 253    | 356    |
| Test 2  | 74     | 281    | 364    |
| Test 3  | 86     | 287    | 372    |
| Test 4  | 80     | 289    | 418    |
| Average | 76     | 278    | 378    |

Table 3-1 Test Results, gallons/Cycle

### Unit 1

Of the three units tested, Unit 1 consumed the least amount of water, averaging 73 percent less than water consumed by Unit  $2^1$  and 80 percent less than water consumed by Unit  $3^2$  (during testing cycles). The reduction in water consumption is mainly attributable to the recycling of water used to create the chamber vacuum. This vacuum is needed for the majority of the conditioning period (9-25 minutes) and for the duration of the drying period (31-32 minutes).

The recirulation pump associated with the Water Mizer system ran for approximately 42 minutes of each test cycle, consuming an average of 0.67 kW or 0.47 kWh. Considering a range of electricity costs between 4 and 15 ¢/kWh this equates to between 1.9 and 7.1 ¢/cycle. The water saved by this recirculation system, compared to Unit 2, was an average of 202 gallons/cycle. Considering water/sewer costs between 4 and 10 \$/kgal, the recirculation system saves between 81 to 202 ¢/cycle. The power consumption of the Water-Mizer system when the pump is not running is limited to controls circuits and was measured at less than 0.1 amps. Therefore, the cost to operate the pump is negligible when compared to the value of water saved with its use.

## Unit 2

Unit 2 still consumed considerably less water during cycles than Unit 3. This is due in large part to the fact that there was no constant drain tempering water flow present in the Unit 2 configuration. The jacket drain of Unit 2 was returned to the boiler plant condensate system, and the chamber drain was tempered only as necessary by the Water-Mizer system.



### Unit 3

Unit 3 was the only test configuration that utilized constant tempering water flow rates for jacket and chamber drains. This tempering water is a significant consumer and a wasteful means of protecting sewer drain temperatures. In fact, Unit 3 consumed more water in standby mode between Tests 3 and 4 (592 gallons), than it consumed during the Test 4 cycle (418 gallons). Furthermore, in a typical operating year for this sterilizer 1,240,900 gallons of water would be consumed in standby mode while 795,400 gallons would be consumed during Dart Test, Leak Test, and Sterilization Cycles<sup>3</sup>. This calculation is based upon a standby water flow rate of 3.1 gpm, 400,285 minutes of standby time, 365 Dart Tests and Leak Tests each, and 1,827 Sterilization Cycles.

The tempering water flow rate, which is manually adjusted by two needle valves near the vacuum venturi (see configuration in Appendix C), was unneccessarily high for the majority of time. During periods of jacket drain only, the discharge temperature was only slightly above the domestic water inlet temperature of 49°F. The non-regulated tempering water flow in this configuration is particularly wasteful during standy periods. If the Water-Mizer system had been installed on Unit 3, the standby water consumption would have been substantially reduced.

### **Other Comments**

The operational cycle water savings benefits of the Water-Mizer (including the recirculation system) illustrated by this testing are two-fold: the reduction of water used to create a vacuum and the reduction of water used to temper discharge flows. The difference in cycle water consumptions between Units 1 and 2 (73 percent) quantify the savings of the first benefit<sup>1</sup>. The difference in cycle water consumptions between Units 1 and 3 (80 percent) quantify the savings of the second benefit as a maximum case (because jacket tempering present in Unit 3 is not neccesary for Unit 1 due to piping differences)<sup>2</sup>. Therefore, the cycle water savings associated with the installation of the Water-Mizer equipment can be concluded to lie somewhere between these two values (73-80 percent).

Based on these test results, if the Water-Mizer (including recirculation) was installed on Unit 3, cycle-related water savings would be anticipated near the high end of the 73 to 80 percent range. Additional savings associated with standby periods would increase overall water savings over the 80 percent mark, possibly exceeding 90 percent<sup>3</sup>. It is also likely that standby savings alone would exceed the cycle savings alone.

The water savings benefit of the installation of the Water-Mizer on Unit 2, was the reduction in chamber discharge tempering. Chamber discharge tempering is only necessary during operational cycles and is therefore minimal when compared to the jacket drain discharge tempering necessary on units that dump jacket condenstate to drain because this flow is necessary through standby periods. Therefore, the benefits of the Water-Mizer are maximized when it is installed on units that dump jacket condensate to drain.



Based on data presented in Table 3-1 the water consumption of this equipment during the Sterilization Cycle can be illustrated as shown in the figure entitled *Characterization of Water Consumption for Sterilizers Tested* contained in Appendix E. In this figure, the blue portion of the bar represents the water flow necessary to temper the sterilizer discharge adequately. This also represents the total cycle flow (80 gallons) if the Water-Mizer with recirculation system is installed (Unit 1). The red portion of the bar represents the quantity of water that is consumed generating a vacuum, and applies if the recirculation system is not installed. The total of the red and blue portions represents the total cycle flow (280 gallons) if the Water-Mizer without the recirculation system is installed (Unit 2). The yellow portion of the bar represents the total of the red, blue, and yellow portions represents the total cycle flow (380 gallons) if the Water-Mizer system is not installed (Unit 3).

As mentioned above and as demonstrated by Unit 3, it is likely that, over prolonged periods, water consumed during standby periods will substantially exceed the water consumed during operational cycles. The greatest potential benefit of the Water-Mizer system, not considering the recirculation system, is the reduction of standby losses in these circumstances.



Conclusions from this testing are as follows:

- The installation of the Water-Mizer including the recirculation system can reduce Sterilizer Cycle water consumption between 73<sup>1</sup> and 80<sup>2</sup> percent. Comparable savings are anticipated for other cycles such as the Dart Test and Leak Test.
- The installation of the Water-Mizer including the recirculation system can reduce total sterilizer water consumption by more than 80 percent. Although not directly tested in this evaluation, these reductions could exceed 90 percent<sup>3</sup>.
- The addition of the recirculation system to the Water-Mizer can reduce Sterilizer Cycle water consumption by 73<sup>1</sup> percent. Comparable savings are anticipated for other cycles such as the Dart Test and Leak Test.
- The installation of the Water-Mizer can drastically reduce standby losses.
- Depending upon the frequency of use, the manual valve settings, and the piping configuration of the condensate drains, sterilizer standby water losses will likely exceed the amount of water consumed during operational cycles (over prolonged periods).
- The cost of operating the electric pump associated with the Water-Mizer recirculation system is negligible when compared to the cost of water saved with its use.
- Considering water and sewer costs to be 5.0 \$/kgal, electricity costs to be 6.0 ¢/kwh, and the total sterilizer water consumption savings to be 90 percent, the installation of the Water-Mizer equipment on Unit 3 at this location would save the owner 9,000 \$/yr in utility costs<sup>4</sup>. This savings represents a typical application of a Water-Mizer installation.



## Appendix A

### The 3 different cycles of the tested sterilizers

| === L E A K T E S T ==<br>CYCLE START AT 11:10:05P<br>ON 4/03/05   | CYCLE START AT 12:02:00A<br>ON 4/04/05   | P R E V A C<br>CYCLE START AT 1:03:07P<br>PM 404/05   |
|--|--|---|
| CYCLE COUNT NL 27933<br>OPERATOR NAC 02<br>STERILIZER VAC 02   | CYCLE COUNT (27934<br>OPERATOR (27934)<br>STERILIZER VAC 02  | CYCLE COUNT<br>OPERATOR<br>STERILIZER VAC 02  |
| U=inHs<br>- TIME T=F P=psis<br>C 11:10:05P 148.2 0P  | STER TEMP = 270.0F<br>CONTROL TEMP = 273.0F<br>STER TIME = 3.5 MIN<br>DRY TIME = 1 MIN   | STER TENN = 278.87<br>CONTROL TEMP = 273.0F<br>STER TIME = 5 MIN<br>DRV TIME = 30 MIN   |
| C 11:11:04P 248.2 19P<br>C 11:12:30P 210.7 15U<br>C 11:12:53P 256.9 26P<br>C 11:14:17P 191 2 27P   | V=inHs<br>- TIME T=F P=psis  | V=inH9<br>- TIME T=F P≃psig   |
| C 11:14:13P 191.2 230<br>L 11:14:51P 270.2 33P<br>L 11:24:50P 165.7 24U<br>L 11:26:51P 163.7 24U<br>L 11:36:50P 158.0 24U<br>LEAK RATE IS:<br>0.0 mmH9/min | C 12:02:00A 151.7 0P<br>C 12:03:00A 260.4 20P<br>C 12:04:26A 210.0 16U<br>C 12:04:45A 249.3 26P<br>C 12:06:06A 194.1 23U<br>C 12:06:25A 241.6 26P<br>C 12:07:44A 194.1 24U | C ~~1:03:07P 131.7 0P<br>C 1:04:07P 245.7 19P<br>C 1:07:59P 220.6 10U<br>C 1:08:46P 258.8 26P<br>C 1:10:57P 213.8 14U<br>C 1:11:35P 251.5 26P<br>C 1:13:37P 214.4 18U |
| L 11:36:51P 158.0 240<br>Z 11:38:10P 165.3 10  | C 12:08:03A 242.2 26P<br>C 12:09:22A 194.2 250<br>S 12:09:55A 270.3 32P  | C 1:14:20P 254.2 26P<br>C 1:16:05P 215.9 190<br>S 1:17:23P 270.1 33P  |
| LOAD 040305<br>TOTAL CYCLE = 0:28:06   | S 12:10:55A 273.7 30P<br>S 12:11:55A 272.8 33P<br>S 12:12:55A 273.2 33P<br>F 12:13:25A 273.0 33P   | \$ 1:19:24P 272.8 34P<br>\$ 1:21:24P 271.9 33P<br>E 1:22:23P 271.6 32P<br>5 1:27:91P 294 5 75   |
|  | E 12:13:46A 231.7 3P<br>E 12:14:45A 186.2 23U<br>Z 12:16:01A 182.2 1V  | E 1:23:01P 224.6 3P<br>E 1:52:59P 182.1 21V<br>Z 1:54:13P 169.3 1V  |
|  | LOAD 040401  | 1000 0404 <u>04</u>   |
|  | TEMP MAX#273.8F<br>TEMP MIN=270.3F   | TEMP MAX=273.1F<br>TEMP MIN=270.1F  |
|  | CONDITION = 7:55<br>STERILIZE = 3:30<br>EXHAUST = 2:36<br>TOTAL CYCLE =14:01   | CONDITION =14:17<br>STERILIZE = 5:00<br>EXHAUST =31:51<br>TOTAL CYCLE =51:08  |
|  | = READY TO UNLOAD =  | = READY TO UNLOAD =   |

# **Appendix B**

|              | Unit      | 1 SJH<br>& 1 Leak ever                  | y 24 hours (appro<br>Water<br>Stand  | age<br>ox. 12:0<br>Meter<br>by Wat | Chart<br>00 a.m.)<br>Installed |              | 3/3/2005<br>2.7 GPM  |
|--------------|-----------|---|--|------------------------------------|--------------------------------|--------------|--|
|              |           | Leak                                    | Test   |                                    | Preve                          | <u>uc #1</u> | and services and   |
| Dart 1       | 22628     | Start                                   | 33787  | S                                  | tart                           | 3397         | 0  |
| start        | 33030     | Finish                                  | 33965  | F                                  | inish                          | 3433         | 6  |
| inish        | 33700     | Total                                   | 178  | Т                                  | otal                           | 36           | 6 , /  |
| l'otal       | 20260     | Cycle CT.                               | 30361  | С                                  | ycle CT.                       | 3036         | 2 3/9/05   |
| Cycle C1.    | 30300     | Cjui en                                 | and the second |                                    |                                |              |  |
| Date         | Wate      | er Meter Gal                            | (  | yele C                             | ount                           |              |  |
| 3/0/05       | 41754     |   | 30372  | in L                               | 1:100                          | IEA.         | DARTILEA   |
| 3-11-05      | 44580     | The second second                       | 30375  | TY .                               | 1:15P                          | ITA.         | DART/LEAD  |
| 3-14-05      | 60937     |   | 30392  | "17                                | 12:000                         | BEA.         | DARTILEAD  |
| 3-16-06      | 67387     |   | 30399  | 7                                  | 2:000                          | IZA          | Destherk   |
| 3-16-05      | 71698     | an ann an Alland                        | 30405  | 6                                  | 8 AM                           | TEA          | DART/LEA   |
| 3-17-65      | 78901     | and the second second                   | 30414  | 9                                  | 1:300                          | 1BA          | DARH/LEA   |
| 3-18-05      | 83725     | The state of the state                  | 30421  | 7                                  | 8 A '                          | IEA          | DARTILEA   |
| 2-21:05      | 100523    |   | 30437  | 16                                 | INOUA                          | 3EA          | OPAT/LEAR  |
| 3-22-05      | 106798    |   | 30445  | 8                                  | 12:00 000                      | 124          | OPET/LENK  |
| 3-23-05      | 113490    |   | 30453  | 8                                  | 3:300                          | 124          | ** *   |
| 3-24-05      | 11950492  | · A A A A A A A A A A A A A A A A A A A | 30461  | 8                                  | 4212                           | 154          |  |
| 325.05       | 123 156   |   | 30466  | 5                                  | The                            | Ira          | at 11  |
| 3-28-03      | 1400730   |   | 30485  | 19                                 | 11:20A                         | 300          | 12 P   |
| 329.05       | 146020    | and a stranger the                      | 30492  | 7                                  | 8:30A                          | 100          |  |
| 3-31-05      | 159527    |   | 30508  | 16                                 | 4:20                           | ZeA          | a 4  |
| 4-1-0        | \$ 163320 |   | 30513  | 5                                  | BAM                            | IEA          | 11   |
| 4-4-05       | 179294    |   | 30529  | 16                                 | 8:38 A                         | 3ca          | X.C. Contraction   |
| 4-4-05       | Recipe.   | YSTEM                                   | 30530  | 1                                  |                                |              |  |
|              | 179359    | <u>access</u>                           | 30531  | 1                                  |                                |              |  |
|              |           |   |  |                                    |                                |              |  |
|              |           | Contraction and the                     |  |                                    |                                |              |  |
| -            |           |   |  |                                    |                                |              |  |
| Carl In Mark | -         |   | -  |                                    |                                |              | the second s |
|              |           |   |  |                                    |                                |              |  |
| A December   |           |   |  |                                    |                                |              |  |
|              |           |   |  |                                    |                                |              |  |
|              |           |   |  |                                    |                                |              |  |

|           | 1 Dart o | & 1 Leak ever<br>Water | y 24 hours (appr<br>r Mizer & Water<br>Stand | ox. 12<br>Meter<br>Iby Wa | :00 a.m.)<br>Installed:<br>ater Usage |              | 2      | 3/8/2005<br>2.8 GPM |
|-----------|----------|------------------------|--|---------------------------|---------------------------------------|--------------|--------|---------------------|
| Dart 7    | est      | Leak                   | Test   |                           | Preve                                 | <u>ac #1</u> |        |                     |
| Start     | 1443.5   | Start                  | 1523   |                           | Start                                 | 15           | 3.5    |                     |
| Finish    | 1523     | Finish                 | 1623   |                           | Finish                                | 40           | 11.5   |                     |
| Total     | 79.5     | Total                  | 100  |                           | Total                                 |              | 248    | bla                 |
| Cycle CT. | 27785    | Cycle CT.              | 27786  |                           | Cycle CT.                             | 27           | 119 :  | 110                 |
| Data      | Water    | Meter Gal              |  | Cycle                     | Count                                 |              |        |                     |
| 240/021   | 100 3890 |                        | 27795  | 9                         | 4:10P                                 | IEA          | LEA    | KIDAR               |
| 21100     | 5524     |                        | 27802  | 7                         | 1.15P                                 | IEA          | LEPA   | 1 ppe               |
| 3-14-05   | 7350     |                        | 27811  | 9                         | 12 '002m                              | 3EA          | LEAR   | loner               |
| 2-15-05   | 8980     |                        | 27818  | 1                         | 2. PM                                 | ILA          | LRAI   | (DART               |
| 2-16-05   | 10334    |                        | 27824  | 6                         | BAM                                   | 1 en         | LE     | K/DAR               |
| 3-17-06   | 12237    |                        | 27832  | 8                         | 1:300                                 | len          | (en)   | -/par               |
| 3-18-05   | 13565    |                        | 27838  | 6                         | 8 A                                   | 199          | LOAL   | DART                |
| 3-21-05   | 16981    |                        | 27854  | 14                        | 11:00 m                               | 3cg          | Local  | 1 par               |
| 3-22-05   | 18363    |                        | 27860  | 6                         | 12 pm                                 | 100          | Courte | 1 Dart              |
| 3-23-05   | 20 4 84  |                        | 27869  | 9                         | 3'30 p                                | JeA          | 8.0    | 41                  |
| 3-24-05   | 22109    |                        | 27876  | 7                         | 4:22                                  | 120          | -      | 1                   |
| 2-25-05   | 22858    |                        | 27880  | 4                         | 70-                                   | 100          | 11     |                     |
| 3-28-05   | 20491    |                        | 27897  | <u>[]</u>                 | 11-10A                                | 3ca          | 40     |                     |
| 2-21-05   | 30452    |                        | 27702  | 6                         | 1.20                                  | 1-1019       |        |                     |
| 4.1.06    | 31498    |                        | 27416  | 17                        | 4120                                  | 2 en         | 44     |                     |
| 4-4-05    | 34476    |                        | 27525  |                           | 6001                                  | 141          |        |                     |
|           |          |                        | e 1 1 2 3                                    | 13                        | C 204                                 | 2004         |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |
| A LANGER  |          |                        |  |                           |                                       |              |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |
|           |          |                        |  |                           |                                       |              |        |                     |

|                  | t<br>1     | Jnit 3 SJH<br>Dart & 1 Leak every | Water U<br>24 hours (app<br>Wat<br>Stan | rox. 12:<br>er Mete<br>dby Wa | cnant<br>:00 a.m.)<br>r Installed<br>iter Usage |                  | 3/3/2005<br>2.8 GPM |
|------------------|------------|-----------------------------------|---|-------------------------------|---|------------------|---------------------|
|                  |            | Leak                              | Test                                    |                               | Preva   | ic #1            |                     |
| Dart To          | <u>est</u> | Start                             | 30381                                   |                               | Start   | 19960            |                     |
| itart            | 30268      | Finish                            | 30551                                   |                               | Finish  | 20338            |                     |
| inish            | 30372      | Total                             | 170                                     |                               | Total   | 378              | abbe                |
| fotal            | 17018      | Cycle CT.                         | 17919                                   |                               | Cycle CT.                                       | 17893            | 5/4/05              |
| Cycle C1.        | 1/910      |                                   |   |                               |   |                  |                     |
| Date             |            | Water Meter Gal                   | and the second second                   | Cycle                         | Count   | Chief            | the -               |
| stolast          | ROU        | 38597                             | 17929                                   | 10 C                          | 1:10 pm   | icr ur           | THINKI .            |
| the ost          |            | 43889                             | 17936                                   | 7                             | 1:15 Pm   | ICA LE           | KI DART             |
| 2-14-05          |            | 54842                             | 17942                                   | 6                             | 12:00 20  | 3EAL             | KAX/DART            |
| 3-15-05          |            | 60923                             | 17948                                   | 6                             | 2 Pm  | I an /           | call parts          |
| 3-16-05          |            | 65478                             | 17954                                   | 6                             | 8 mm  | 1 en li          | EAK/PORT            |
| 3-17-05          |            | 72710                             | 17962                                   | 8                             | 1:30 p  | I CA L           | AR DART             |
| 3-18-05          |            | 77356                             | 17968                                   | 6                             | 87  | 1=+ 4            | ear / Padt          |
| 3-21-05          |            | 73200                             | 17975                                   |                               | Incom   | SCA LA           | at locks            |
| 3-22-05          |            | 98921                             | 17981                                   | 6                             | 12 pm   | 100 b            | enk / User          |
| 3.23.05          |            | 105592                            | 17988                                   | 7                             | 3:302   | ICA.             |                     |
| 3.24.05          |            | 111883                            | 1976                                    | <u> </u>                      | 4 42  | 125              | 11 24               |
| 3-25-05          |            | 115 607                           | 18001                                   | 10 Da                         | 1 days  | 2                | 33                  |
| 3-20-05          | -          | 13/11/2                           | 18012                                   | - 6                           | 8:30  | Jeh              |                     |
| 3.27.03          | +          | 1490.94                           | 18026                                   | 13                            | 4:20  | 2 - 0            | 25 25               |
| 4-1-00           | 4          | 152434                            | 18029                                   | 2                             | 8.4   | 100              |                     |
| 4-4-06           | -          | 16202.1                           | 18035                                   | 6                             | S 38 A  | 340              |                     |
|                  |            |                                   | 1                                       |                               |   |                  | West Marshall       |
|                  |            |                                   |   |                               |   |                  |                     |
|                  |            |                                   |   |                               |   |                  |                     |
|                  |            |                                   | Contraction of the second               |                               |   |                  |                     |
|                  |            |                                   | -                                       | and the second second         |   |                  |                     |
|                  |            |                                   |   |                               |   |                  |                     |
|                  |            |                                   |   |                               |   |                  |                     |
|                  |            |                                   |   |                               | Sector Sector                                   |                  |                     |
| The second       |            |                                   |   |                               |   | Mar and a second |                     |
| A State of State |            |                                   |   |                               |   |                  |                     |
|                  |            |                                   |   |                               |   |                  |                     |

## **Appendix C**







## **Appendix D**

CYCLE START AT 3:86:23P 0N 4/05/05 V=inH9 Pepsig 188.1 89 2333.4 179 2593.5 189 2593.5 189 2593.5 189 2588.8 269 268.8 269 268.8 269 2796.1 211 2786.1 319 273.8 379 273.8 379 273.8 379 273.8 379 273.8 379 273.8 19 272.8 379 272.8 10 \* = READY TO UNLOAD = 040507 = 278.0F = 273.0F = 5 MIN = 38 MIN CONDITION = 0:25:29 STERILIZE = 0:05:00 EXMAUST = 0:32:14 TUTAL CYCLE = 1:02:43 38546 UAC B1 Т=F TEMP MAX=273.2F TEMP MIN=270.1F STER TENP = CONTROL TENP = STER TIME = DRY TIME = CYCLE COUNT OPERATOR \_\_\_\_\_ STERILIZER **JULE** LUAD . 1 IU UNLOAD = 17:33:12P 17:33:12P 47:05705 V≕in49 P=psia 111111 ļ READY TO UNLOAD 040506 VAC 01 30539 270.0F 273.0F 5 min 30 min CONDITION =19:52 STERILIZE = 5:00 EXHAUST =32:02 TOTAL CYCLE =56:54 148.3 239.8 2296.6 2219.1 2219.1 2219.1 2209.8 2273.1 2273.1 2273.1 2273.1 2273.1 156.7 1176.7 1176.7 TEMP MAX=273.2F TEMP MIN=270.1F Ï 11 U B B <del>8</del>8 STER TEMP : CONTROL TEMP : STER TIME : DRY TIME : - TIME - TIME - TIME - TIME - 12:33:12P - 12:35:13P - 12:55:13P - 12:55:12P OPERATOR \_\_\_\_ STERILIZER START CYCLE COUNT CYCLE LOAD ===== 10: 55:350 4/05/05 U=inHs P=⊳sia ZEADY TO UNLOAD ≠ 105.9 0P 2338.2 19P 219.2 10U 2558.8 26P 2558.8 26P 261.9 26P 261.9 26P 261.7 14U 261.7 14U 273.1 32P 273.1 32P 273.1 32P 273.1 32P 273.1 32P 273.1 22U 273.1 22P 273. 040505 = 279.0F = 273.0F = 5 MIN = 30 MIN 39538 UAC 01 =18:37 = 5:00 =32:07 LE =55:44 TEMP MAX=273.2F TEMP MIN=270.1F 1=F ₽40 STER TEMP = CONTROL TEMP = STER TIME = GRY TIME = CONCITION = STERILIZE = EXHGUST = TOTAL CYCLE = CVCLE COUNT OPERATOR \_\_\_\_\_ STERILIZER START 10:35:359 10:35:349 10:44:96 10:45:969 10:45:996 10:45:996 10:45:289 10:45:179 10:59:129 10:59:119 10:55:119 10:55:119 11:89:949 11:59:924 11:59:924 11:59:169 TIME CVCLE LOAD I. .===== 9:19:13A 4/05/05 V=i∩Hs P=Psis 040504 276.0F 273.0F 5 min 30 min 30537 91 153.0 153.0 2646.8 2612.9 2612.9 273.0 275 UAC = 9:07 = 5:00 =31:40 =45:47 TEMP MAX=273.2F TEMP MIN=270.4F 1 ппц 19 19 11 STER TEMP = CONTROL TEMP = STER TIME = DRY TIME = 9:19:13A 9:20:13A 9:20:13A 9:22:51A 9:22:51A 9:22:51A 9:24:45A 9:26:10A 9:26:10A 9:25:20A 9:27:20A 9:33:20A 9:30A 9:33:20A 9:34:20A 10:34:20A 10:34:20A 10:34:20A 10:34:20A 10:34:20A 10:34:20A 10:34:20A CYCLE COUNT OPERATOR \_\_\_\_\_ STERILIZER CONDITION = STERILIZE = EXHAUST = TOTAL CYCLE = START ------TIME CYCLE LOAD 

#### Unit #1 Test Log Printouts

3:40:05P 4/05/95 V≖ínH9 P=⊭sí9 ===== b E C d C ===== 840567 - 279.9F - 273.9F 5 Min 38 Min 27947 62 141.8 241.6 2259.1 2259.1 2259.1 2251.9 2251.9 2272.4 2272.4 2272.4 2272.4 2272.4 2272.4 2772.1 2772.1 2772.4 2777.4 27777.4 277 CONDITION =13:19 TERILIZE = 5:80 EXHAUST =31:49 FOTAL CVCLE =59:09 READY TO UNLO' 3 TEMP MAX=272.6F TEMP MIN=270.1F Ĩ 11 81 14 14 YCLE START AT TIME YCLE COUNT PERATOR \_\_\_\_\_ STERILIZER STER ONTROL STER DRY TIME 000 1:36:06P 4/05/05 V=inH∋ P=⊳sig VDV TO UNLOAD = 17200 1720 1720 1720 1720 17200 1720 17200 17200 17200 17200 940506 270.0F 273.0F 5 Min 30 Min 27946 18 =13:19 = 5:00 =31:50 =50:89 UAC  $\odot$ MAX=273.1F MIN=270.1F L L 0.8.0.0 START AT ON STER TEMP = CONTROL TEMP = STER TIME = DRV TIME = ====== P R E CYCLE COUNT OPERATOR \_\_\_\_ STERILIZER CVCLE 1:35:05P 1:37:05P 1:42:29P 1:42:14P 1:42:14P 1:45:35P 1:45:35P 1:45:35P 1:45:35P 1:45:35P 1:45:35P 1:45:35P 1:55:25P 2:25P 1:55:25P 2:25P 1:55:25P 2:55:25P CONDITION STERILIZE EXHAUST TOTAL CVCL TEMP TIME CVCLE LOAD ı. Ċ V=inHs P≍psis EVCLE START AT 18:16:07A ===== P R E V A C ===== H 040505 27945 273.0F 273.0F 5 MIN 30 MIN UNLOAD 140.9 242.7 216.1 216.1 216.1 219.5 =12:56 = 5:00 =31:55 =49:51 TEMP MAX=273.0F TEMP MIN=270.1F ц. П 1 SONDITION = TERILIZE = EXHAUST = OTAL CYCLE = 2 19:16:079 19:17:379 19:17:379 19:29:559 19:25:594 19:25:594 19:25:594 19:25:594 19:25:594 19:25:594 19:25:594 19:25:594 19:25:594 19:25:594 19:25:594 19:25:564 19:34:034 19:35:564 11:85:564 READY CONTROL 1 CONTROL 1 STER 1 DRV 1 TIME CLOAD CLOCKED COLLEGE COLLEGE A CONTRACTOR OF 9:08:58A 4/05/05 V=inHs P=₽sist = 270.0F = 273.0F = 5 MIN = 30 MIN 27944 UAC 02 READY TO UNLOAD <u>640564</u> 136.1 248.3 248.3 2609.3 1260.7 2260.7 2243.8 2243.8 2243.8 2243.8 2243.8 2243.8 2272.6 100.7 2772.6 100.7 2772.6 100.7 2772.6 100.7 2772.6 2772.6 2772.7 2777.7 2775.7  $\odot$ = 8:15 = 5:00 =31:37 =44:52 æ • MAX=273.4F • MIN=270.1F Ξ. TEMP TIME TIME CYCLE COUNT OPERATOR \_\_\_\_\_ STERILIZER 9:08:58A 9:09:57A 9:11:526 9:11:526 9:11:526 9:11:526 9:11:526 9:12:156 9:15:156 9:15:156 9:15:156 9:15:136 9:15:136 9:15:137 9:15:137 9:15:137 9:52:54 9:52:524 CVCLE STAK, CONDITION = STERILIZE = EXHAUST = TOTAL CYCLE = STER -CONTROL 1 STER 1 DRY 1 3011 TEMP LOAD T I 

#### Unit # 2 Test Log Printouts

| <br>                     | CYCLE COUNT 18047<br>OPERATOR                             | STER TEMP = 270.0F<br>CONTROL TEMP = 273.0F<br>STER TIME = 5 MIN<br>DRV TIME = 30 MIN | utinHe<br>- Time Tef Pepsis | C 2:51:14P 117.8 BP<br>C 2:52:13P 236.8 17P<br>C 2:57:12P 210.4 10U<br>C 2:58:05P 258.4 26P                             | C 3:80:14P 211.7 11U<br>C 3:80:54P 254.0 26P<br>C 3:02:51P 209.1 15U<br>C 3:02:51P 209.1 15U   | C 3:85:337 239.4 267<br>C 3:85:337 239.4 267<br>S 3:86:339 278.2 337<br>S 3:18:399 273.3 337<br>S 3:18:399 273.2 337<br>E 3:11:399 273.4 337<br>E 3:12:319 223.2 37<br>E 3:42:281 172.5 230<br>F 3:42:281 172.5 230    | L0AD 840507                         | TEMP MAX≈273.6F<br>TEMP MIN≈270.2F | CONDITION =15:26<br>STERLLIZE = 5:00<br>EXHAUST =5:11<br>TOTAL CYCLE =52:37   |
|--------------------------|---|---|-----------------------------|---|--|--|-------------------------------------|------------------------------------|---|
| CYCLE START AT 10:55:01A | CVCLE COUNT 18046<br>OPERATOR UAC 03<br>STERILIZER UAC 03 | STER TEMP = 270.0F<br>CONTROL TEMP = 273.0F<br>STER TIME = 5 MIN<br>DRY TIME = 30 MIN | U=irH9<br>- TIME T=F P=p≤is | C 10:55:01A 127.7 1P<br>C 10:57:01A 127.2 1P<br>C 10:57:02A 209.0 1EU<br>C 10:59:32A 209.0 1EU<br>C 10:59:32A 209.0 1EU | C 11:01:130 202.0 170<br>C 11:01:41A 248.2 26<br>C 11:03:14A 248.2 26<br>C 11:03:14A 262.0 190 | U 11:05:2474 250.0 269<br>0 11:05:228 209 200<br>3 11:05:228 201 34<br>3 11:08:136 273.4 34<br>5 11:10:138 273.4 34<br>5 11:11:129 273.7 33<br>E 11:11:129 273.7 33<br>E 11:11:1496 275.7 33<br>7 11:41:466 1781 240   | LOAD 840505                         | TEMP MAX=273.7F<br>Temp Min=270.1F | CONDITION =10:12<br>STERILIZE = 5:00<br>EXHAUST =31:56<br>TOTAL CYCLE =47:00  |
| CVCLE START AT 9:45:17A  | CYELE COUNT 19045<br>OPERATOR<br>STERILIZER UAC 03        | STER TEMP = 270.0F<br>CONTROL TEMP = 273.0F<br>STER TIME = 5 MIN<br>DRY TIME = 30 MIN | y=inH3<br>- T[ME T=F P=psig | C 9:45:174 126.2 1P<br>C 9:45:16A 126.2 1P<br>C 9:46:16A 247.3 21P<br>C 9:47:56A 205.4 12U<br>C 9:48:76A 257 2 26P      | C 9:50:204 197.4 190<br>C 9:50:27A 249.1 26P<br>C 9:51:58A 200.9 200                           | C 9:52:24A 249.0 26P<br>C 9:53:55A 205.9 21U<br>S 9:55:44A 270.1 34P<br>S 9:58:44A 273.6 34P<br>S 9:58:44A 273.6 34P<br>S 9:59:43A 273.6 34P<br>E 9:59:43A 273.1 34P<br>E 10:00:15A 138.8 23U<br>E 10:00:15A 138.8 23U | z 19:31:32H 142.3 10<br>LOAD 040505 | TEMP MAX=273.9F<br>TEMP MIN=270.1F | CONGITION = 9:27<br>STERILIZE = 5:00<br>EXHAUST =5:1:51<br>TOTAL CVCLE =46:18 |
|                          | SYCLE COUNT 18944<br>DPERATOR 18944<br>STERILIZER 0AC 03  | STER TEMP = 270.0F<br>CONTROL TEMP = 273.0F<br>STER TIME = 5 MIN<br>DRV TIME = 30 MIN | U=inHa<br>- TIME T=F P=psig | C 8:51:32A 133.6 0P<br>C 8:52:32A 248.9 21P<br>S 8:54:02A 188.2 19U<br>S 8:54:230 188.2 19U                             | 8:55:490 181.8 220<br>8:55:496 181.8 220<br>8:57:326 185.4 220<br>9:57:326 185.4 220           | 8:55:528 248.2 26P<br>8:59:158 187.2 220<br>8:59:554 270.1 34P<br>9:03:554 270.1 34P<br>9:03:548 270.1 34P<br>9:04:548 273.1 34P<br>9:05:288 273.1 33P<br>9:05:288 220.6 37P<br>9:05:288 125.1 23U                     | . 2.20.244 144.4 14                 | TEMP MAX=274.0F<br>TEMP MIN=270.1F | CONDITION = 8:22<br>STERILIZE = 5:00<br>CAHAUST = 31:42<br>'DTAL CYCLE =45:04 |

### Unit # 3 Test Log Printouts

# **Appendix E**











# Appendix F



Unit 1





Unit 3



### The Sterilizers Tested



Water Meter (Typical of 3)



**Power Consumption Measurement** 



Water-Mizer Recirculation Tank Installed on Unit 1

## Appendix G

| Footnote<br>Calculation | Cycle Water Savings Low Range Result                                 |
|-------------------------|--|
| 1                       | Average Water Consumption Unit 1, per cycle                          |
|                         | Average Water Consumption Unit 2, per cycle<br>278 gallons           |
|                         | Average Water Savings with Retrofit (Unit 1 vs Unit 2)<br>73 percent |

| Footnote    | Cycle Water Savings High Range Result                  |  |  |
|-------------|--|--|--|
| Calculation |  |  |  |
| 2           | Average Water Consumption Unit 1, per cycle            |  |  |
|             | 76 gallons   |  |  |
|             | Average Water Consumption Unit 3, per cycle            |  |  |
|             | 378 gallons  |  |  |
|             | Average Water Savings with Retrofit (Unit 1 vs Unit 3) |  |  |
|             | 80 percent   |  |  |

#### Measured Parameters and Assumptions

| Assumed Water/Sewer Rate                            | 5.00 \$/kgal        | Measured kW draw of pump   | 0.67 kW       |  |
|---|---------------------|--|---------------|--|
| Assumed Electricity Utility Rate                    | 0.06 \$/kWh         |  |               |  |
| Typical Dart Test Water Consumption, Pre-Retrofit   |                     | Typical Dart Test Cycle Time                                       |               |  |
| 113 gallons   |                     | 30 minutes   |               |  |
| Typical Leak Test Water Consumption, Pre-Retrofit   |                     | Typical Leak Test Cycle Time                                       |               |  |
| 174 gallons   |                     | 13 minutes   |               |  |
| Typical Sterilization Cycle Water Consumption, Pre- | Retrofit            | Typical Sterililzation Cycle Time                                  |               |  |
| 378 gallons   |                     | 60 minutes   |               |  |
| Assumed Number of Dart Tests per Year               |                     | Entering Domestic Water Temperature                                |               |  |
| 365 (once per day)                                  |                     | 49 degrees F   | 17.05 Btu/lb  |  |
| Assumed Number of Leak Tests per Year               |                     | Drain Discharge Temperature during Standby Pre-Retrofit            |               |  |
| 365 (once per day)                                  |                     | 54 degrees F   | 22.05 Btu/lb  |  |
| Assumed Number of Sterilization Cycles per Year     |                     | Drain Discharge Temperature during Standby Post-Retrofit           |               |  |
| 1827 (7 per workday/261 workdays pe                 | r year)             | 140 degrees F  | 107.98 Btu/lb |  |
| Cycle Water Savings with Retrofit (maximum of calc  | ulations one above) | Jacket Condensate Drain (Saturated Liquid at atmospheric pressure) |               |  |
| 80 percent  |                     | 212 degrees F  | 180.81 Btu/lb |  |
| Typical Standby Water Flow Pre-Retrofit             |                     |  |               |  |
| 3.1 average gpm                                     | 1,562 average lb/h  |  |               |  |

| Footnote    | Annual Water Consumptions and Estimated Savings  |  |  |  |  |
|-------------|--|--|--|--|--|
| Calculation |  |  |  |  |  |
| 3           | Calculated Time Spent in Cycles (Total Number of Dart, Leak, and Sterilization Cycles multiplied by the time for each) 125,315 minutes   |  |  |  |  |
|             | Calculated Standby Time per year (Total Minutes per year, 525,600 minus times spent in cycles)<br>400,285 minutes  |  |  |  |  |
|             | Calculated Flow of Jacket Condensate Drain During Standby (before tempering) (via heat balance)<br>49 average lb/h 0.097 average gpm   |  |  |  |  |
|             | Calculated Standby Water Flow Post-Retrofit (via heat balance)<br>39 average lb/h 0.077 average gpm  |  |  |  |  |
|             | Calculated Annual Water Consumption, Pre-Retrofit (Dart, Leak, and Sterilization Cycles plus Standby Losses)<br>795,361 gallons during cycles<br>1,240,884 gallons during standby<br>2,036,245 gallons total |  |  |  |  |
|             | Calculated Annual Water Consumption, Post-Retrofit (Dart, Leak, and Sterilization Cycles plus Standby Losses)<br>189,894 gallons   |  |  |  |  |
|             | Calculated Annual Water Savings Resulting from Retrofit<br>1,846,351 gallons<br>1,846 kgal<br>91 percent<br>0.230 delore   |  |  |  |  |
|             | 5,230 U0idis   |  |  |  |  |

| Footnote    | Annual Utility Savings Considering Cost to Operate Pump  |
|-------------|--|
| Calculation |  |
| 4           | Pump Operating Time For Each Cycle (the pump was observed to operate 42 minutes of each 60 minute cycle)<br>70 percent |
|             | Calculated Annual Pump Operating Time<br>1,462 hours   |
|             | Calculated Annual Pump Power Consumption<br>980 kWh<br>59 dollars  |
|             | Calculated Annual Overall Utility Savings (Water Utility Savings minus Costs to Operate Pump)<br>9,171 dollars         |